DATASHEETQuickFix Elastic

1. Identification of the substances/preparation and of the company/undertaking

Date issued 11.10.2004
Revision 06.05.2011
Product name Quick Fix Elastic
Index no. Version 1.2

Use of the substance/preparation Plaster QuickFix Elastic

Producer

Company name Plum A/S

Postal address Frederik Plums Vej 2

Postcode 5610
Place name Assens
Countryl Danmark
Tel. 0045 6471 2112

 Fax
 0045 6471 21 25

 E-mail
 info@plum.dk

 Website
 http://www.plum.dk

2. Backing

Composition	CD elastic fabric (65% cotton, 35% polyamide)
Colour	Pink
Weight (g/m ²)	80-100

>30

3. Adhesive

01710110	
Composition	Hypoallergenic acrylic polymer
Weight (g/m2)	55-65

4. Coated backing

MD Tensile strenght (N/cm)

Weight (g/m ²)	135-165
Adhesion to steel by peeling at 180° (N/cm)	>1,2
M.V.T.R. (g/m ² /24h/37°C)	>500
Water resistant	no

5. Absorbent pad

Total weight (g/m²)	100-160
Absorbency (g/100 cm ²)	9-16

Absorbent layer Viscose and polypropylene nonwoven

Colour White

Antiadherent layer Microperforated PE film

6. Liner

Monosiliconated paper

Weight (g/m ²)	57-63
Colour	White

QuickFix Elastic

7. Primary packing

Individual pouch: paper + clear film

Paper composition Cold sealed paper

Paper weight (g/m^2) 40-50

Film composition Polyethylene terephthalate / polyethylene

Film weight (g/m2) 52-67

8. Handling, storage and instructions for use

Further information To be used before expiry date.

Single-use.

QuickFix Elastic is a medical devices according to Council Directive

93/42/EEC of 14 June 1993 concerning medical devices.

No instructions for use are needed for this device as it can be used

safely without such instructions

9. Document changes

06-05-2011 KIS Section 8 and 9 inserted